

Patient Instructions

An AFM™ Ag dressing has been applied to your wound(s). The dressing contains an antimicrobial ingredient to help reduce the risk of infection. It may be covered with additional dressings to absorb excess drainage and provide added protection.

The dressing has been applied in the appropriate direction for your wound:

- To remove excess fluids with the printed side up
- To add or retain moisture, with the printed side down

It is recommended that you:

- Leave all dressing materials undisturbed until your next appointment
- Change the secondary gauze dressings as needed, leaving the AFM Ag dressing in place.

Notes: _____

If the dressing becomes soiled or damaged, contact our office.

Clinician: _____

Date: _____